



Kinetic Systems, Inc.

Quote Request for ELpF

Please print this form, fill in desired selections, and fax to Kinetic Systems Inc. at **617.522.6323**

ELpF Lightweight Benchtop Vibration Isolation System

Select desired accessories below.

Faraday Cage

Ordering Information

ELpF -

Capacity

- 01 = 100 lb Capacity
- 02 = 300 lb Capacity

Top Description

- 11 = White plastic laminate
- 12 = Ferromagnetic stainless steel
- 13 = Black plastic laminate
- 14 = White antistatic laminate
- 22 = Ferromagnetic stainless steel with 1/2"-20 x 1" centers (SPILLPRUF)
- 25 = Ferromagnetic stainless steel with 1/2"-20 x 0.5" centers (SPILLPRUF) MaxxAlign
- 32 = Ferromagnetic stainless steel with M6 holes x 25mm centers (SPILLPRUF)
- 35 = Ferromagnetic stainless steel with M6 holes x 12.5mm centers (SPILLPRUF) MaxxAlign

*Custom tabletop sizes available

PURCHASE ORDER NUMBER _____

DESCRIPTION OF EQUIPMENT TO BE ISOLATED _____

EQUIPMENT SIZE AND WEIGHT _____

ESTIMATED CENTER OF GRAVITY _____

NAME _____ TITLE _____

COMPANY _____

ADDRESS _____ CITY _____ STATE _____ ZIP _____

TEL: _____ FAX: _____

E-MAIL ADDRESS _____