



Kinetic Systems, Inc.

Quote Request for Series 9100

Please print this form, fill in desired selections, and fax to Kinetic Systems, Inc. at **617.522.6323**

Series 9100 Workstation Ordering Information

91 - -

Gross Load Capacity (Including Tabletop) 01 = Up to 800 lbs (364kg) 02 = Up to 1300 lbs (591kg)	Tabletop Description	Weight lbs/sq.in (kg/sq.m)
	00 = No Tabletop	
	21 = 2" (51mm) Composite w/white plastic laminate	.14* (98)
	22 = 2" (51mm) Composite w/stainless laminate	
	23 = 2" (51mm) Composite w/black plastic laminate	
	24 = 2" (51mm) Composite w/antistatic laminate	
	30 = 4" (101mm) VT Gray granite	.40 (281)
	40 = 2" (51mm) Honeycomb - no holes	.12 (84)
	45 = 2" (51mm) SPILLPRUF Honeycomb 1/4-20 x1" centers	
	46 = 2" (51mm) SPILLPRUF Honeycomb M6 x 25mm centers	
	80 = 4" (101mm) Honeycomb - no holes	.15 (105)
	85 = 4" (101mm) SPILLPRUF Honeycomb 1/4-20 x1" centers	
	86 = 4" (101mm) SPILLPRUF Honeycomb M6 x 25mm centers	

* 30" x 72", 36" x 60", and 36" x 72" composite tops weigh .11 lb/sq. in.

Select desired KSI accessories below.

Suggested accessories:

- GR-L Guard Rails Front & Rear (L = Table Length)
- PA2 Padded Armrests (pair) *
- RC-4 Retractable Casters (set of 4)

Additional accessories:

- CL10 Class Ten Cleanroom Compatible Finish
- EOS Electrical Outlet Strip
- ENC Enclosure**
- FC Faraday Cage**
- FS Fixed Shelf (24" or 30" Deep Tables)
- FSC Fixed Shelf Custom (tables deeper than 30")
- GLU Gross Load Upgrade (800lbs. to 1300lbs)
- HE-5 Height Extension (5" Integral)
- HER-5 Height Extension (5" Field Retrofit)
- HOP Hole Pattern (Stainless Steel Laminate Only)
- KS Keyboard Shelf* (Fixed: 10" x 24")
- KF Kinematic Feet (set)
- MS Monitor Support (Mounts to any Corner)
- DMS Monitor Support (Double)
- ORES Overhead Rear Equipment Shelf
- AERS Overhead Rear Additional Rear Shelf
- OSES Overhead Side Equipment Shelf
- AESS Overhead Side Additional Side Shelf
- PAC Portable Air Compressor
- SR Seismic Restraint
- SDR Side Rails
- SS-10 Sliding Shelf* - 10" Wide
- SS-12 Sliding Shelf* - 12" Wide
- SS-14 Sliding Shelf* - 14" Wide
- SSC Sliding Shelf*, Custom (>14" Wide, > 30" Deep)

* Requires Guard Rails

**Requires Guard Rails and Side Rails

PURCHASE ORDER NUMBER _____

DESCRIPTION OF EQUIPMENT TO BE ISOLATED _____

EQUIPMENT SIZE AND WEIGHT _____

ESTIMATED CENTER OF GRAVITY _____

NAME _____ **TITLE** _____

COMPANY _____

ADDRESS _____ **CITY** _____ **STATE** _____ **ZIP** _____

TEL: _____ **FAX:** _____

E-MAIL ADDRESS _____