



Kinetic Systems, Inc.

Quote Request for Series 8002

Please print this form, fill in desired selections, and fax to Kinetic Systems, Inc. at **617.522.6323**

Series 8002 Active Isolation Table

Ordering Information

8002 - -

Tabletop Size

Depth **Width**

1 = 24" (609mm) 1 = 36" (914mm)
 0 = 30" (762mm) 2 = 48" (1219mm)
 2 = 36" (914mm) 3 = 60" (1524mm)
 3 = 48" (1219mm) 4 = 72" (1828mm)
 4 = 60" (1524mm)

Gross Load Capacity (only 2 Capacities)
 600 lbs (273kg)
 1200 lbs (545kg)

Tabletop Description

40 = 2" (51mm) Honeycomb/no holes
 45 = 2" (51mm) SPILLPRUF Honeycomb 1/4-20 x 1" centers
 46 = 2" (51mm) SPILLPRUF Honeycomb M6 x 25mm centers
 47 = 2" (51mm) SPILLPRUF Honeycomb
 MaxxAlign 1/4-20 x 0.5" centers
 48 = 2" (51mm) SPILLPRUF Honeycomb
 MaxxAlign M6 x 12.5mm centers
 80 = 4" (101mm) Honeycomb/no holes
 85 = 4" (101mm) SPILLPRUF Honeycomb 1/4-20 x 1" centers
 86 = 4" (101mm) SPILLPRUF Honeycomb M6 x 25mm centers
 87 = 4" (101mm) SPILLPRUF Honeycomb
 MaxxAlign 1/4-20 x 0.5" centers
 88 = 4" (101mm) SPILLPRUF Honeycomb
 MaxxAlign M6 x 12.5mm centers

Select desired accessories below.

Suggested accessories:

- GR-L Guard Rails Front & Rear (L = Table Length)
- PA2 Padded Armrests (pair) *
- RC-4 Retractable Casters (set)

Additional accessories:

- CL10 Class Ten Cleanroom Compatible Finish
- EOS Electrical Outlet Strip
- ENC Enclosure**
- FC Faraday Cage**
- FS Fixed Shelf (24" or 30" Deep Tables)
- FSC Fixed Shelf Custom (tables deeper than 30")
- HOP Hole Pattern (Stainless Steel Laminate Only)
- KS Keyboard Shelf* (Fixed: 10" x 24")
- MS Monitor Support (Mounts to any Corner)
- DMS Monitor Support (Double)
- ORES Overhead Rear Equipment Shelf
- AERS Overhead Rear Additional Rear Shelf
- OSES Overhead Side Equipment Shelf
- AESS Overhead Side Additional Side Shelf
- SR Seismic Restraint
- SS-10 Sliding Shelf* - 10" Wide
- SS-12 Sliding Shelf* - 12" Wide
- SS-14 Sliding Shelf* - 14" Wide
- SSC Sliding Shelf*, Custom (>14" Wide, > 30" Deep)

* Requires Guard Rails
 ** Requires Guard Rails and Side Rails

Note: Active Isolators Are Available Without Workstations. Consult your KSI Sales Representative.

PURCHASE ORDER NUMBER _____

DESCRIPTION OF EQUIPMENT TO BE ISOLATED _____

EQUIPMENT SIZE AND WEIGHT _____

ESTIMATED CENTER OF GRAVITY _____

NAME _____ **TITLE** _____

COMPANY _____

ADDRESS _____ **CITY** _____ **STATE** _____ **ZIP** _____

TEL: _____ **FAX:** _____

E-MAIL ADDRESS _____